

KINEO MULTISTATION

TECHNICAL SPECIFICATIONS

Dimensions (LxWxH)

260x90x185 cm

Footprint Area (LxWxH)

305x90x185 cm

Weigth

332 Kg

Electric motor

Nominal power 1.6kW

Power supply

Tension: 230V - Frequency: 50/60Hz

Maximal current: 8A

Maximum cable extraction

4 m

Min/Max Leg Extension load

10-630 N/m

Min/Max Pulley load

1-110 Kg

Min/Max Squat load

1-220 Kg (up to 440 Kg with dedicated accessory)

Min/Max Isokinetic Speeds

0.1-4 m/s

Min/Max Leg Extension Isokinetic Speeds

12-778 °/s

FEATURES

- Wi-Fi connection and remote control
- 2 15.6" touch-screen displays
- 2 USB ports
- Integrated memory and automatic archive system
- Italian high-quality leather
- CE 0476 certification

HARDWARE

- Electric motor that produces variable and controlled torque
- Pulley workstation
- Leg Extension workstation
- Squat workstation

ACCESSORIES





- Kit Leg Curl
- Kit Leg Press
- Lat Machine
- Squat wooden platform

MEDICAL COMPLIANCE



KINEO MULTISTATION

SOFTWARE CONFIGURATIONS

| TRAINING METHODS | VERSION 3.0 | VERSION 5.0 | VERSION 7.0 |
|---|-------------|-------------|-------------|
|  ISOTONIC NIT (NO INERTIA) | ✓ | ✓ | ✓ |
|  ISOTONIC WIT (WITH INERTIA) | ✓ | ✓ | ✓ |
|  ELASTIC | | ✓ | ✓ |
|  SMART VISCOUS | | ✓ | ✓ |
|  VARIABLE LOAD | | | ✓ |
|  ISOMETRIC | | | ✓ |
|  ISOKINETIC | | | ✓ |

| TESTS | VERSION 3.0 | VERSION 5.0 | VERSION 7.0 |
|--|-------------|-------------|-------------|
| CONCENTRIC MAX STRENGTH TEST | ✓ | ✓ | ✓ |
| CONCENTRIC BALANCE STRENGTH TEST | ✓ | ✓ | ✓ |
| POWER TEST | ✓ | ✓ | ✓ |
| ECCENTRIC MAX STRENGTH TEST | | ✓ | ✓ |
| ECCENTRIC BALANCE STRENGTH TEST | | ✓ | ✓ |
| JUMP POWER TEST | | ✓ | ✓ |
| ISOMETRIC MAX STRENGTH TEST | | ✓ | ✓ |
| ISOMETRIC BALANCE STRENGTH TEST | | ✓ | ✓ |
| ISOKINETIC MAX STRENGTH TEST | | | ✓ |
| ISOKINETIC BALANCE STRENGTH TEST | | | ✓ |
| ISOKINETIC CONCENTRIC/ECCENTRIC RATIO TEST | | | ✓ |